

Special-Use Bags



Capacity is measured in fluid ounces when the bag is closed and tab is folded over three times. Volume and dimensions are approximate; bags should not be used at temperatures above 180° F (82° C). Bags can be frozen to any temperature, even liquid nitrogen — careful handling required after freezing.

Thio-Bags® and Coli-Test Bag for Water Sampling



BO1040WA
With sodium thiosulfate.
100 ml
3" W x 7 1/4" L
(7.5 x 18.5 cm)
2.25 mils
(.057 mm) thick
Box of 100

BO1254WA
With sodium thiosulfate.
300 ml
4 1/2" W x 9" L
(11.5 x 23 cm)
2.5 mils
(.064 mm) thick
Box of 100

BO1440WA
30 mg sodium thiosulfate, meets NELAC Standard.
100 ml
3" W x 7 1/4" L
(7.5 x 18.5 cm)
2.25 mils
(.057 mm) thick
Box of 100

Stand-Up Bags



BO1402WA
With sodium thiosulfate.
100 ml Stand-Up
3" W x 7 1/4" L
(7.5 x 18.5 cm)
2.5 mils
(.064 mm) thick
Box of 100

BO1403WA
With sodium thiosulfate.
300 ml Stand-Up
4 1/2" W x 9" L
(11.5 x 23 cm)
3.0 mils
(.076 mm) thick
Box of 100

BO1417WA
Coli-Test bag for color-changing reagents.
100 ml Stand-Up
3" W x 7 1/4" L
(7.5 x 18.5 cm)
2.5 mils (.064 mm) thick
Box of 100

Homogenizer Blender Bags



BO1018WA
13 oz. (384 ml)
5 1/2" W x 7 1/4" L
(13 x 19 cm)
3.0 mils
(.076 mm) thick
Box of 500

BO1063WA
24 oz. (710 ml)
6" W x 9" L
(15 x 23 cm)
4.0 mils
(.102 mm) thick
Box of 500

BO1196WA
24 oz. (710 ml)
6" W x 9" L
(15 x 23 cm)
4.0 mils
(.102 mm) thick
Box of 500

BO1195WA
55 oz. (1,627 ml)
7 1/2" W x 12" L
(19 x 30 cm)
4.0 mils
(.102 mm) thick
Box of 500

BO1239WA
52 oz. (1,538 ml)
7" W x 13" L
(18.5 x 33 cm)
4.0 mils
(.102 mm) thick
Round bottom seal
Box of 500

Spoon/Scoop Bags



BO1478WA
With sterile teaspoon.
18 fl. oz. (532 ml)
4 1/2" W x 9" L
(11.5 x 23 cm)
2.5 mils
(.064 mm) thick
Box of 100

BO1350WA
With 2-oz. (60 ml) scoop.
18 fl. oz. (532 ml)
4 1/2" W x 9" L
(11.5 x 23 cm)
2.5 mils
(.064 mm) thick
Box of 50

Bubble Bags



BO1483WA
Inside dimensions:
2 1/4" (5.7 cm) dia.
4" W x 8" L
(10.2 x 20.3 cm)
Pack of 10

BO1484WA
Inside dimensions:
6" (15.2 cm) dia.
9" W x 12" L
(30.5 x 22.9 cm)
Pack of 10

Nuclear Bags



BO1407WA*
Holds dosimeter, extra long tabs.
3 1/4" W x 7" L
(9.5 x 18 cm)
2.25 mils
(.057 mm) thick
Box of 500

BO1485WA*
Holds dosimeter, extra long tabs.
5" W x 10" L
(12.5 x 25.4 cm)
2.5 mils
(.064 mm) thick
Box of 500

NEW!
BO1520WA*
Holds dosimeter, extra long wireless tabs.
5" W x 10" L
(12.5 x 25.3 cm)
2.5 mils
(.064 mm) thick
Box of 500

NEW!
BO1521WA*
Holds dosimeter, extra long wireless tabs.
6" W x 10" L
(15 x 25.4 cm)
2.5 mils
(.064 mm) thick
Box of 500